

# Tranquility COLLECTION




shown in  
Hazelnut

Copper



Stocked with heat  
and massage

Hazelnut



Sterling



Calypso



Shiraz



Silt



Mahogany



The UC 480 is going hybrid! This model will continue to offer a traditional biscuit-style backrest and plush chaise pad, but is now sized to fit a broader height range of individuals. This medium/large chair can now be everyone's favorite seat in the house.

## UC 480MLA

*UltraComfort*<sup>™</sup>  
AMERICA  
Positioning you for life.

PROUDLY  
  
MADE IN AMERICA  
with Domestic & Imported Materials

# Tranquility COLLECTION

## UC 480

Model Number	UC480
Size	Medium/Large
Weight Capacity	375 LB
Recline Positions	3
Chaise Seat	Yes
Pocketed Coil Seating System	No
Back Type	Biscuit
Head Rest Cover	Yes
Arm Covers	Yes
Battery Back Up	Yes
Overall Width	35"
Overall Height	41.5"
Width Between Arms	22"
Floor to Top of Seat	22"
Seat Depth	22"
Seat to Top of Back	31"
Seat to Top of Arm	7"
Seat Angle Raised	20°
Distance Required from Wall Reclined	19"
Extension Reclined from Back of Seat to End of Foot Rest	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Frame, Recline Mechanism	Yes
Warranty on Electrical parts	3 YR
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	3 YR
Warranty on Electrical Labor	1 YR
Warranty on Fabric and Foam, and Pocketed Coils	1 YR
Removable Back Option	Yes
Customizable *	No
Left Hand Control Option Available	Yes
Extra Pocket Option Available	Yes
4-Motor Heat & Massage Available	No
2-Motor Heat & Massage Available	Yes
Foot Rest Extension Option Available	Yes
Smart Tek Standard	Yes

\* See Customize Chart for dimension limits  
 Sizes may vary depending on fabric, filling material, upholstery and carpet thickness.  
 All measurements taken on concrete floor using levels and metal rulers.



### User-friendly feature of every UltraComfort lift and recline chair:

Easy to use battery backup system. Each UltraComfort lift and recliner has a battery backup system in the electrical transformer which operates on two 9 volt batteries. In the event of a power failure, the battery system will raise the chair from the reclined to the seated position. The batteries should be replaced after one use.



UltraComfort™  
 AMERICA

Positioning you for life.

UltraComfort  
 Phone: 866-793-5100  
 Fax: 866-255-9962  
 web: [www.ultracomfort.net](http://www.ultracomfort.net)



Literature is current at time of printing. UltraComfort reserves the right to make changes to the product or literature at any time.



Complies with UL 73 and CSA C22.2 No. 68

Listing # E112417  
 FNAE 80-214A

