

Slip into the casual comfort of the all-new Montage Collection! Our new design features a plush, generously stuffed seam back, better seat support and a more casual armrest design for hours of relaxation. In four different sizes and six standard colors, the Montage collection is sure to become a household favorite!



Montage 542-L For people 5'11" -6'2"





Montage UC 542-L

| Model Number | UC542-L |
|---|---------|
| Size | Large |
| Weight Capacity | 375 LB |
| Recline Positions | 3 |
| Chaise Seat | Yes |
| Pocketed Coil Seating System | No |
| Back Type | Seam |
| Head Rest Cover | Yes |
| Arm Covers | Yes |
| Battery Back Up | Yes |
| Overall Width | 35.5" |
| Overall Height | 46" |
| Width Between Arms | 23" |
| Floor to Top of Seat | 21" |
| Seat Depth | 21" |
| Seat to Top of Back | 30" |
| Seat to Top of Arm | 8" |
| Seat Angle Raised | 15° |
| Distance Required from Wall Reclined | 20" |
| Extension Reclined from Back of Seat to End of Foot Rest | 44.5" |
| Life Time Warranty, Chair Frame, Lift Frame, Recline Frame, Recline Mechanism | Yes |
| Warranty on Electrical parts | 3 YR |
| Pro Rated Years 4 Thru 7 | Yes |
| Warranty on Mechanical Labor | 3 YR |
| Warranty on Electrical Labor | 1 YR |
| Warranty on Fabric and Foam, and Pocketed Coils | 1 YR |
| Removable Back Option | No |
| Customizable * | Yes |
| Left Hand Control Option Available | Yes |
| Extra Pocket Option Available | Yes |
| 4-Motor Heat & Massage Available | Yes |
| 2-Motor Heat & Massage Available | No |
| Foot Rest Extension Option Available | Yes |
| Smart Tek System | Yes |

^{*} See Customize Chart for dimension limits Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



User-friendly feature of every UltraComfort lift and recline chair:

Each UltraComfort lift and recliner has a battery backup system in the electrical transformer which operates on two 9 volt batteries. In the event of a power failure, the battery system will raise the chair from the reclined to the seated position. The batteries should be replaced after one use.





Positioning you for life.

UltraComfort Phone: 866-793-5100 Fax: 866-255-9962 web: www.ultracomfort.net



